

## Increased Viability and Genomic Integrity of CRISPR-modified hiPS cells selected with WOLF<sup>®</sup> Cell Sorter Microfluidic Technology


Elizabeth Roberts<sup>1</sup>, Allison Songstad<sup>1</sup>, Elena Rubio de la Torre<sup>2</sup>, Emily Bozek<sup>2</sup>, Archana Shankar<sup>1</sup>, Gabriela Goldberg<sup>1</sup>, Florian Krach<sup>1</sup>, Stefan Aigner<sup>1</sup>, Gene Yeo<sup>1</sup>, Lawrence Goldstein<sup>1</sup>. 1. Sanford Consortium for Regenerative Medicine, UCSD. 2. NanoCelect Biomedical, Inc., San Diego, CA

**Background:** Human iPS cells (hiPSC) can be modified utilizing CRISPR/Cas9 gene editing technology to generate *in-vitro* models of human diseases. Critical shortcomings in the generation of CRISPR-modified clonal lines are: 1) low viability/proliferation, and 2) genomic abnormalities following selection via fluorescent activated cell sorting. To address these issues, we compared microfluidic cell sorting (WOLF<sup>®</sup> Cell Sorter, NanoCelect) to a traditional electrostatic droplet based cell sorter in Alzheimer's Disease (AD) and Autism Spectrum Disorder (ASD) model cell lines. We evaluated colony formation and karyotypes after sorting.


**Methods:** Human fibroblasts were reprogrammed into hiPSC using retroviruses containing the four Yamanaka Factors. Resulting hiPSCs were then CRISPR-modified to knock-down AD- or ASD-related genes. GFP was used as reporter for selection. Cells were sorted in parallel with the WOLF<sup>®</sup> and a traditional droplet-based cell sorter. After sorting, cells were plated on MEFs and grown for 1-2 weeks, picked, expanded, and frozen. Colonies that recovered after cryogenic storage were submitted for microarray karyotype analysis.

**Results:** ASD-related SHANK3 and AD-related SWE, V717F, and ICE cell lines displayed improvements in viability and genomic integrity following microfluidic sorting on the WOLF<sup>®</sup>, compared to a traditional droplet cell sorter. The number of proliferating colonies, normal karyotypes, as well as overall freeze-thaw viability were all increased on the WOLF<sup>®</sup> microfluidic platform compared to a droplet sorter.

### Autism Spectrum Disorder hiPCs



### Alzheimer's Disease hiPCs



**Figure 1: Viability and genomic integrity after sort:** **A** Number of colonies that proliferated after sorting and plating for the SHANK3-knockdown CRISPR-modified cells. Numbers above bars indicate number of cells sorted and plated. **B** SHANK3 karyotyped colonies that displayed genomic abnormalities. **C** Number of colonies that proliferated after sorting for the SWE, V717F, and ICE -positive CRISPR modified cells. Numbers inside bars indicate number of cells sorted and plated (x1000). **D** Number of colonies that recovered after cryogenic freeze-thaw measured via trypan blue staining. **E** Number of SWE, V717F, and ICE karyotyped colonies that displayed genomic abnormalities.